

Bringing Personalized Medicine to the Clinic



cMET: Simultaneous and Quantitative Analysis of Copy Number Variation and Gene Expression in a Single Reaction

Molecular Med Tri-Con
Feb., 13, 2013

PrimeradX

Lilly Kong, DVM
CSO, PrimeradX

PrimeradX – Simplifying Molecular Diagnostics with a Proprietary Platform

Only technology available to enable Multi-Modal Multiplex testing

- Ability to detect and quantify DNA, RNA, mRNA and miRNA all in a single-well reaction
- Proprietary technology with large IP estate
- Automated platform that dramatically improves lab workflow, reduce costs, and enables new markets

Technology advantages provide opportunity for differentiated content; several assays in development

- cMET/EGFR and other oncology assays in partnership with Pharma
- KRAS/NRAS/BRAF, EGFR and several ID products currently in development

Technology is validated and in the market now

- Eli Lilly companion diagnostics relationship, with multiple other opportunities in progress
- Initial open platform placements ramping up

Technology has unique and large commercial opportunity

- “Open platform” approach
- Cleared IVD kit sales
- Companion diagnostic partnerships

Regulatory path is clear and in process



Unique Technology Addresses Current Problems While Enabling New Markets and Applications


The Marriage of PCR and Capillary Electrophoresis

Integrated PCR, CE, fluidics, detection and analysis software

Real-time detection of PCR products separates targets by size

Sampling as the reaction progresses allows quantitation

Multiplex, multimodal, quantitative results provide unique, clinically actionable data




The ICEPlex System – Fully Automated, Real-Time, Multiplex qPCR

Walk-away Workflow with Automated Reporting of Assay Results

Thermal
Cycler



On-board Reagents




Capillary
Cartridge

- Dual lasers, multiple dye sets, dozens of targets per dye
- Assay dynamic range (and simultaneous detection) of 10 – 10,000,000 copies of multiple targets in a single sample
- Innovative software to track, analyze and report results
- Proven reliability – customer experience
- Manufactured under QSR, ISO and GMP standards
- Flexible software: User-definable assay conditions for LDT capabilities, company-developed assay design software speeds assay/product development

There is No Comparable Technology Available

ICEPlex Enables Real-time Multi-modal, Multiplex, Quantitative Tests




Next Generation
Clinical Tests

Multi-Modal Testing is the Future of MDx – Breakthrough for Labs and Clinicians

Diagnostics Now Require Next Generation “All in one Well” Capabilities

- End-point Qualitative Detection
- Real-time Quantification
- SNP
- Fusion Gene
- Gene copy #
- mRNA
- microRNA
- Controls



A real-time, multiplex, multi-modal “All in one well” MDx solution

Only PrimeraDx Can Do This

Unique, Critical Capability in the Molecular Dx Space

Many Tests, One Reaction/One Platform

Current practice: many platforms, many techs



PrimeraDx



ICEPlex

Oncology Multi-Modal Panel:

Fusion gene variants

+

SNPs

+

Gene Expression Signature


+

MicroRNAs

+


Standards and controls

29 Targets plus 9 Standards and Controls



96 Patient Samples

6 days



1 morning

Regulatory Plan and Quality Systems Enable PDx in the Market

Market Needs Drive Approach

Achieving Regulatory Clearances


- IVD kit sales are a major portion of market
- CDx Partners need a Cleared Platform
 - Will be a major differentiator for PrimeraDx



Regulatory Milestones

- Operate under QSRs ✓
- ISO 13485 ✓
- CE marks ✓
- FDA *de novo* 510(k) 2013


PrimeraDx Progress on Quality Systems and Regulatory Approvals



PrimerDx's Open Platform Product: Next Generation qPCR

Users Design Complex Multiplex Assays with Ease....

In Silico Multiplex Assay Design SW Tools Enable Customers



Automated Workflow Simplifies Testing Procedures

.....and PrimerDx Provides the Tools Needed to Diagnose Disease

PrimeraDx Open Platform Products

Instrument

ICEPlex enables CLIA labs to run laboratory developed tests *and* PrimeraDx IVDs



Consumables



Four cartridge sizes
8 and 24 well
48 well
96 well

Universal Reagent Kit
Capillary Electrophoresis plates
All on-board consumables

Products

Open Platform Product

Instrument, software and consumables needed for CLIA labs to run LDTs



PrimeraDx Technology Has Very Broad Clinical IVD Utility

Proven Capabilities in All Major Diagnostic Areas....

Oncology

- Expression panels
- Copy number variation
- SNP panels
- Polysomy
- Insertions
- Deletions
- Fusion products (replaces FISH)
- Methylation
- Combinations of **any and all** of the above

More in Development



Others

- Food Safety Testing
- Pharma QC
- Drug Metab
- Genetic Dis

More in Development

ID

- Quantitative, multi-pathogen detection
- Mixed pathogen panels
- Viral load
- Resistance detection
- Hospital acquired infections
- Multiple sample/swab types
- Panels that represent physician ordering patterns

More in Development

Panels for demonstration purposes only. Not for clinical diagnostic use.

PrimeraDx Products Have Very Broad Clinical IVD Utility

Proven Capabilities in All Major Diagnostic Areas....

Oncology

- KRAS/NRAS/BRAF
- EGFR
- cMET/EGFR CNV with cMET Expression
- cMET Mutation Panel
- EML4-ALK
- Lymphoma
- BCR-ABL
- Methylation
- microRNA/mRNA/gDNA Panel

More in Development



Others

- Food Safety Testing
- Pharma QC
- Drug Metab
- Genetic Dis

More in Development


ID

- Fungal Panel
(direct blood detection)
- Transplant Panel
(multiple viral loads)
- C. difficile*
- STI Panel
(yeast, virus, bacteria, parasites)
- Ultra-Sensitive JCV/BKV
(differentiation, viral loads)
- Respiratory Panel
(pathogens including drug-resistant strains)




More in Development

Panels for demonstration purposes only. Not for clinical diagnostic use.

c-MET/EGFR Copy
Number Variation
and cMET Gene
Expression
Single-Tube Assay



c-MET and EGFR – Critical in Many Pathways, Targeted by Many Drugs



Guo A et al. PNAS 2008;105:692-697

Why c-MET???

cMET is a receptor tyrosine kinase that, after binding with its ligand, hepatocyte growth factor (HGF), activates many signaling pathways, driving proliferation, motility, migration and invasion.

Although c-MET is important in the control of tissue homeostasis under normal physiological conditions, it has also been found to be aberrantly activated in human cancers via mutation, amplification or protein overexpression. Dysregulation and constitutive activation of c-MET leads to cell proliferation, cell survival, angiogenesis, invasion and metastasis.


Overexpression/amplification of c-MET has been observed in various carcinomas, including gastric, NSCLC, colorectal, kidney tumors.

Amplification of c-MET has been identified as one mechanism to confer resistance to EGFR-specific tyrosine kinase inhibitors in lung cancers.

Challenges of Current Platforms & Multi-modality

IHC: Protein Expression/phosphorylation

- ❑ Low sensitivity of antibodies
- ❑ strong background staining, weak target antigen staining and autofluorescence



FISH: Gene/Chromosomal CNV

- ❑ Inter-lab discordance 20%+
- ❑ 5-7 day turn-around time



Real-Time PCR: Gene Expression / SNP/CNV

- ❑ Multiplex limitations on most other PCR platforms

Many biomarkers in many complex pathways targeted by many drugs:

- ❑ Creates a need for single tube, multiplex and multi-modal tests



A “Killer App” For a High Unmet Need – Only PrimeraDx Can Do This

cMET & EGFR COPY NUMBER VARIATION

gDNA cMET, EGFR and Ref Gene #1 on Chromosome 7



gDNA Ref Gene #2



gDNA Ref Gene #3



Primers located in introns selective for gDNA

ONE TUBE ASSAY



- Multi-modal
- Multiplex (18+6 plex)
- Quantitative
- FFPE
- Whole process controls

QUANTITATIVE cMET GENE EXPRESSION PROFILE

mRNA cMET and Ref Gene #1 on Chromosome 7



mRNA Ref Gene #2



mRNA Ref Gene #3





Primers spanning exon/exon junctions selective for mRNA

Detects and Measures

- CNV
- Gene Expression
- Polysomy
- Reference genes


Panels are designed for demonstration purposes only. Not for clinical diagnostic use.

Single Tube CNV and Gene Expression Analysis – Gastric




Single Tube CNV and Gene Expression Analysis – Gastric

cMET and EGFR CNV




cMET Gene Expression




Single Tube CNV and Gene Expression Analysis – Gastric

Sample	Tissue Source	cMET Copy #	cMET Expression	Ch 7 Polysomy
SNU-1	Gastric	2	No	No
SNU-5	Gastric	>10	High	Yes


cMET and EGFR CNV



cMET Gene Expression




cMET Gene Expression and CNV Analysis



- Data were normalized to reference genes Ref2 and Ref3 and fold-change calculated
- No cMET amplification or expression in SNU-1
- SNU-5 shows significant cMET amplification and overexpression

Panels are designed for demonstration purposes only. Not for clinical diagnostic use.


Chromosome 7 Polysomy Status



No Change

Panels are designed for demonstration purposes only. Not for clinical diagnostic use.

cMET 13-Target Single-tube Mutation Panel



Panels for demonstration purposes only. Not for clinical diagnostic use.

Summary


- An 18-target, 24-plex, single-tube multimodal assay designed to detect amplification of cMET and EGFR genes, expression of cMET and polysomy of chromosome 7, as well as a 13-target single-tube cMET mutation panel have been developed on the ICEPlex system
- Quality result: Built-in controls and redundancy assure quality and precision
- Quicker TAT: Impacts patient care
- All-in-one assay: Simplifies lab operation and improves lab economics
- Minimum nucleic acid input: Resolves specimen size issue
- Versatile sample types: FFPE blocks, Fresh-frozen tissues, Cell lines, Blood

PDx Makes “Liquid Biopsy” Possible – Massive Potential Market in Cancer Screening

KRAS/BRAF Panel: Detects DNA in Serum/Plasma, and Also Works in FFPE

One-tube KRAS / BRAF Panel

- KRAS G12R
- KRAS G12C
- KRAS G12S
- KRAS G12D
- KRAS G12A
- KRAS G12V
- KRAS G13R
- KRAS G13C
- KRAS G13S
- KRAS G13D
- KRAS G13A
- KRAS G13V
- BRAF V600E/D




Screening and disease recurrence monitoring require repeated testing

1% selectivity - highly sensitive serum/plasma assay

Panels for demonstration purposes only. Not for clinical diagnostic use.


Multiplex Detection of EGFR Mutations in One Single Reaction



Panels for demonstration purposes only. Not for clinical diagnostic use.

EML4-ALK: Fusion Gene Assay Detects Eight Fusion Variants

Fusion Detection Will Displace FISH in the Marketplace – Huge Potential




<u>Target</u>	<u>Amplicon size</u>
Var 1	161
Var 2	166
Var 3a	128
Var 3b	154
Var 4a	143
Var 4b	124
Var 5a	133
Var 5b	118
ALK	173
EML4 wt	185

Fusion Gene Variant Detection on the ICEPlex is simple, and provides definitive answers in just a few hours.

Panels for demonstration purposes only. Not for clinical diagnostic use.

BCR-ABL Fusion Gene Assay Detects Nine Fusion Variants



PrimeradX Panel Represents the Best Test Available



Individual plasmid target templates were detected by the single tube BCR-ABL Assay

Panels for demonstration purposes only. Not for clinical diagnostic use.


Lymphoma Gene Expression Panel Differentiates GCB from ABC



- A novel quantitative 19-plex mRNA expression profiling assay designed to allow DLBCL tumor classification on FFPE specimens in a single tube PCR reaction was developed on the ICEPlex system.
- The ICEPlex DLBCL assay allowed discrimination of ABC and GCB cell lines based on specific target expression patterns generated from < 1ng of RNA from FFPE-isolated material.

Panels for demonstration purposes only. Not for clinical diagnostic use.

Simultaneous Quantification of mRNA and miRNA expression



Currently being investigated for complex molecular assays (multimodal multiplex qPCR) supporting drug development with expectation that taking a single assay through regulatory approval will be more feasible

Platform selection is important:

Robust assay performance

Ease of development

Regulatory path

Clinical laboratory accessibility

Meeting the need of the CDx effort

CDx Partnership for IVD Development

Multiyear, Multiproduct Deal with Eli Lilly; Others in Process



Features	<ul style="list-style-type: none"> • Assay design • Assay development • Under Design Control • Tissue-specific 	<ul style="list-style-type: none"> • Final robustness testing • Transfer to manufacturing • cGMP production lots • IDE status • CE Mark 	<ul style="list-style-type: none"> • Trial support to CLIA CROs • V&V studies • Concordance analyses • PMA prep & submission • XUS registrations 	<ul style="list-style-type: none"> • Pre-launch conversion from CROs to clinical labs • Facilitate lab partner reimbursement • Support post-approval surveillance trial
Deliverables	IVD-track assay for use in CLIA or R&D setting	IVD-track IUO for use in registrational trial	<ul style="list-style-type: none"> • PMA submission • Global registrations 	Global distribution


Partner funds product development, expands pipeline of “killer-apps.”

PrimeraDx retains all commercial rights.

PrimerADx – Providing Solutions for Today's Dx Needs

- The only **Fully Quantitative, High Multiplex** qPCR Platform
- A completely **Automated Workflow**
- Allows for **Multi-modal** assay configurations
- The ideal platform for **Companion Diagnostics**
- Numerous complex products in **Diverse Clinical Areas**

Primer Design Use Case: BCR-ABL Panel



ICEPlex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

Assay Design Workflow – Automated Wizard-like User Interface

- Target 1 → Automated Primer Selection
- Target 2 → Automated Primer Selection
- Target 3 → Automated Primer Selection
- Target 4 → Automated Primer Selection
- Target 5 → Automated Primer Selection
- Target 6 → Automated Primer Selection
- Target 7 → Automated Primer Selection
- Target 8 → Automated Primer Selection
- Target 9 → Automated Primer Selection
- Target 10 → Automated Primer Selection



ICEplex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

Primer Design for Target Region

Match all of the following:

- Any Field contains SC5314 mitochondrial
- Organism contains candida albicans

Throttle sequence list (records):

Name	Description	Organism	Taxonomy	Common Name	Accession
NC_002653	Candida albicans SC5314 mitochondrial genome	Candida albicans SC5314	Fungi	mitochondrion (in Candida albicans SC5314 NC_002653)	NC_002653
NW_139443	Candida albicans SC5314 chromosome 4 Cig19-D1158...	Candida albicans SC5314	Eukaryota	Fungi	NW_139443

Sequence View: Doublet (Self) Annotations Text View History Notes

Extract R.C. Translate Allow Editing Add/Edit Annotations Synthesize & Predict Save

Primer Probe dialog box:

- Exclude primer pairs with a single site
- Exclude primers with GC-rich ends
- Maximum number of G or C at 3' end: 6



Design New Primers

Select Task: Design New Design with Existing

Primer design uses Primer3. Please cite [Primer3](#) if you publish results.

Forward Primer DNA Probe Reverse Primer

Region Input Options:

- Included Region: From 1 To 40,420
- Target Region: From 4,660 To 4,830
- Product Size Between: 150 And 250
- Optimal Product Size: 1
- Number of Pairs to Generate: 3

Advanced Options:

- Allow degeneracy: 1
- Calculate dimerization (slow)
- Inverse PCR

Primer

Size Min:	18	Optimal:	20	Max:	27
Tm Min:	57	Optimal:	60	Max:	63
%GC Min:	20	Optimal:	50	Max:	80
Product Tm Min:	0	Optimal:	0	Max:	0

Max Tm Difference: 100 GC Clamp: 0
Max Hairpin Score: 8 Max Primer-Dimer Score: 3
Max Poly-X: 5 Max 3' Stability: 9

Allow primers inside target with penalty: 0


Primer Picking Weights Tm Calculation Settings



ICEplex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

Automated Primer Filtering & ePCR

In silico prediction of off-target amplifications



ICEplex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

The Finalized BCR-ABL Panel – PCR Steps

Simple and Painless High Multiplex Assay Design

PCR Reaction				PCR Primer	uM in 25x n	Sequence
				e13	5	GCTGACCAACTCGTGTGTGAACTCC
	Stock	[Final]		e19	5	TCATGGAGGAGGTGGGCATCTACC
FastStart buffer	10x	1x	5	e6	5	ATATTGCTGTCTCAGGCCAATGCTCAGTTTGC
MgCl₂	25mM	0.5 mM	1	e8	5	CCACCCTGACCACCCTTGCTG
dNTPs	10mM	0.3 mM	1.5	e1	5	TGTCCGAGGCCACCATCGTG
FastStartTaq	5 U/ul	5 U	1	ABLwt	3.75	ACATCACGCCAGTCAACAGTCTG
Calibrators mix	25x	0.3x	0.6	ABL _a 2	7.5	/56-FAM/ACCAACGAGCGGCTTCACTCAGA
Primer mix	30x	1x	2	ABL _a 3	2.5	/56-FAM/ATAAAAATTAATTTTGGTTTGGGCTTCACACCATTCC
Template (cDNA)			14	ABLwt	3.75	/56-FAM/TCTCGGAGGAGACGTAGAGCTTG
H₂O			24.9			
		Total	50			
PCR Protocol						
Initial Denaturation	1x	96C	600 sec			
Pre-CE Cycling	17x	64C	40 sec			
		72C	40 sec			
		96C	10 sec			
CE Cycling	24x	64C	40 sec			
		72C	190 sec			
		96C	10 sec			
CE Sampling	11x	11 separations				
		18, 20...38				
Separation Time		350 sec				
Total Run Time		3hrs 29min				

ICEplex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

The Finalized BCR-ABL Panel – Corresponding Sizes for Each Targets

Simple and Painless High Multiplex Assay Design

Target	Amp Size	Forward					Reverse		
e13a2	128	e13	GCTGACCAACTCGTGTGTGAACTCC				ABLa2	/56-FAM/ACCAACGAGCGGCTTCACTCAGA	
e13a3	136	e13	GCTGACCAACTCGTGTGTGAACTCC				ABLa3	/56-FAM/ATAAAAATTAATTTGGTTTGGGCTTACACCATTCC	
e19a2	140	e19	TCATGGAGGAGGTGGGCATCTACC				ABLa2	/56-FAM/ACCAACGAGCGGCTTCACTCAGA	
e6a2	168	e6	ATATTGCTGTCAGGCCAATGCTCAGTTTGC				ABLa2	/56-FAM/ACCAACGAGCGGCTTCACTCAGA	
e8a2	179	e8	CCACCCTGACCACCCCTTGCTG				ABLa2	/56-FAM/ACCAACGAGCGGCTTCACTCAGA	
e14a2	204	e13	GCTGACCAACTCGTGTGTGAACTCC				ABLa2	/56-FAM/ACCAACGAGCGGCTTCACTCAGA	
e14a3	212	e13	GCTGACCAACTCGTGTGTGAACTCC				ABLa3	/56-FAM/ATAAAAATTAATTTGGTTTGGGCTTACACCATTCC	
e1a2	144	e1	TGTCCGAGGCCACCATCGTG				ABLa2	/56-FAM/ACCAACGAGCGGCTTCACTCAGA	
e1a3	153	e1	TGTCCGAGGCCACCATCGTG				ABLa3	/56-FAM/ATAAAAATTAATTTGGTTTGGGCTTACACCATTCC	
ABLwt	223	ABLwt	ACATCACGCCAGTCAACAGTCTG				ABLwt	/56-FAM/TCTCGGAGGAGACGTAGAGCTTG	

ICEPlex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

Insert the PCR Plate into the Thermal Block




ICEplex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

Insert the PCR Plate into the Thermal Block




ICEPlex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

Input the Platemap and Run Protocol



ICEPlex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

Results!




ICEplex system and assays have not been approved by the FDA for IVD. This information is for demonstration purpose only

BCR-ABL Fusion Gene Assay Detects Nine Fusion Variants

PrimeradX Panel Represents the Best Test Available

One-tube
BCR-ABL SNP
Panel



Individual plasmid target templates were detected by the single tube BCR-ABL Assay

Panels for demonstration purposes only. Not for clinical diagnostic use.

PrimeraDx

The Multiplex PCR Company