

Automated High Multiplex qPCR Platform for Rapid Detection of Fungal Pathogens

Kiran Madanahally Divakar, Jason Lei, Smriti Gupta and Lilly Kong. PrimeraDx, Mansfield, MA 02048

Abstract

Quantitative PCR using probe-based real-time detection is restricted to few targets within a single reaction. The ICEPlex[®] System was developed to offer an automated, high multiplexing solution. The ICEPlex System combines PCR thermal cycling with dynamic, sequential amplicon separation and multi-color quantitative detection by capillary electrophoresis (CE) in a single integrated system. Unlike probe-based qPCR, the ICEPlex System directly measures amplicon accumulation through incorporation of labeled primers.

Here we demonstrate that the ICEPlex System is capable of simultaneous detection and discrimination of fungal pathogens, *Candida albicans*, *C. glabrata*, *C. krusei*, *C. parapsilosis*, *C. tropicalis*, and *Cryptococcus neoformans*, in a single reaction. In addition, the ICEPlex Fungal Assay consists of an extraction control, a sensitivity control, an internal control, a pan-fungal primer set, and a calibration control. Analytical studies indicate 10-20 fungal cell sensitivity per milliliter of whole blood, good dynamic range, and high specificity.

The ICEPlex System is an automated bench top instrument that has overcome many of the limitations of traditional real-time PCR-based multiplexing. Using this 16-plex ICEPlex Fungal Assay (12 set for yeast detection, one set for pan-fungal detection, one set for extraction control detection, one set for calibrator control, and one optional set for internal control detection), we demonstrate the performance of the ICEPlex System in detection of medically important fungal pathogens.




Figure 2. Principles of the multiplex real-time PCR detection on the ICEPlex System

Technology Overview

Unlike any other detection system, ICEPlex integrates a PCR thermal cycler that accommodates a standard 96-well PCR plate, a capillary electrophoresis (CE) system with a replaceable CE cartridge, and fluidic pumps with on-board reagents. The ICEPlex System has two solid state lasers (488 nm and 639 nm excitation) and a spectrophotometer with a CCD camera, allowing two different dyes to be detected simultaneously in single well. (Figure 1)




Figure 1. Components of ICEPlex fully automated real-time PCR and CE system

Unlike other real-time PCR technologies, STAR (Scalable Target Analysis Routine) technology detects and quantifies fluorescently labeled PCR products with unique sizes on the ICEPlex System by sequential sampling and separation on a capillary gel electrophoresis (CE) in real time. Multiple quality controls and calibrators can be embedded in the PCR reaction to ensure quality of the reported results. Built-in software automatically generates amplification curves by converting the fluorescent signals corresponding to the detected PCR amplicons once the run on ICEPlex is complete. A cycle threshold (C_t) for all PCR targets is reported in the Clinical Report (Figure 2).

Materials and Methods

Primer design: A unique region of the genome for each target in the assay was identified by bioinformatics analysis followed by homology search with Basic Local Alignment Search Tool (BLAST) from the National Center for Biotechnology Information (Bethesda, MD). Geneious[™] Pro software (Auckland, New Zealand) was used to carry out, (1) ClustalW analyses (European Bioinformatics Institute, Cambridge, UK) of unique sequences; (2) primer design on the consensus sequences (Steve Rozen and Helen J. Skaletsky (2000) Primer3 on the WWW for general users and for biologist programmers); and (3) primer-primer interaction analyses among different primers (Cross-Hyb, Primera Dx, Inc., Mansfield, MA).

Genomic DNA templates: Genomic DNAs of *Candida*, *Cryptococcus*, *Aspergillus*, and other microbial DNA used in the study were obtained from American Type Culture Collection (Manassas, VA) whereas human genomic DNA was obtained from Roche Applied Science (Indianapolis, IN). Sensitivity and Internal controls were custom synthesized by Integrated DNA Technologies (IDT, Inc., Coralville, IA). Calibrator control is a product from PrimeraDx, Inc (Mansfield, MA).

PCR setup and amplification conditions: The PCR reactions contained: 1X Qiagen[®] Multiplex PCR Master Mix with 1 U of HotStarTaq[®] DNA Polymerase (Qiagen, Germantown, MD), 0.2 μM of each gene-specific forward and reverse primer pairs (IDT, Inc), of which one is labeled with FAM-dye (Table 1), and 0.35 X of the Universal Calibrators (PrimeraDx). Twenty five microliter multiplex PCR reactions were subjected to thermo cycling in a standard 96-well PCR plate on the ICEPlex System. PCR amplification conditions were as follow:

- 98°C for 10 minutes
- 16 cycles at 60°C for 20 sec., 72°C for 30 sec. and 96°C for 20 sec.
- 20 cycles at 64°C for 110 sec., 72°C for 150 sec., 96°C for 10 sec.

Nucleic acid purification: Nucleic acids of yeast cells after lysis were purified/extracted using NucliSENS[®] easyMag[®] according to the manufacturer's instructions.

Primer Name	Primer Sequence	Primer Name	Primer Sequence
>ALBICANS_M_F	ACTGTAGCTGAGCGTAAGACATTAGGTTAT	>CA_RPS_F	ACGTTGATGTTGCCATTC
>ALBICANS_M_R	ATAGTAACCTCCGGGTCCGAAAGGTATAA	>CA_RPS_R	ATATGCTCTTCGACGGCATTTTG
>GLABRATA_R_F	ATATCGTTGGGGAAACGCTCT	>CG_CIT_F	ATCTGCCAAGTCGACCCCTAACG
>GLABRATA_R_R	AGCTTGGCGCTCGTGTGCC	>CG_CIT_R	ATATCCGACCAAGTGTGTGGTGTG
>KRUISIL_M_F	AGCTTTACCCCTCTGAATTTCTG	>CK_ABC_F	ATGGATAGCGAGGACGGCGTG
>KRUISIL_M_R	CCACCGGACGCGTTGT	>CK_ABC_R	GCCTGGTGTACCCATCCGT
>PARAPSILOSI_M_F	TCTGAAGGTTGTACGAAATGGGAAAAA	>CP_SCR3_F	AGACATTGCCAACAAAAGCTCC
>PARAPSILOSI_M_R	CAGAAGACCGTAGTATCGCTGACCAACATTT	>CP_SCR3_R	ATCTGACCAAGCTTGAGCGT
>TROPICALIS_M_F	ATATCAACCATCTCCAGCCCCATC	>CP_CIT_F	AGCCAGCCAAAGCTGAAGAAG
>TROPICALIS_M_R	ATATTTGCTCACCGTGTACCTCCA	>CP_CIT_R	GGGCTTTTGGCAATTCCTTT
>NEOFORMANS_R_F	ATCAACGGATCTCTGGCTCCACAT	>CN_CIT_F	ATGAGGACGCGCTGACGAGAA
>NEOFORMANS_R_R	ACCCAAATCCAAGTCCCAACAGGTAATAAAA	>CN_CIT_R	ATAGCAGTGTATCGGGGGTCT
Pan-fungal ITS34_F	CGTGAATCATCGAATCTTTGAACGCAC	Extraction control_pombe_F	ATGGGTATCCGATTCGCTGGTT
Pan-fungal ITS34_R	GGTTTCTTTTCTCCGCTATTGATATGC	Extraction control_pombe_R	GGGCAACGGTCAATGAGTTC
Internal control_F	GGAAGGCTGTGCATCCAAAGATCCT		
Internal control_R	CATGAATGGTGTGAGCTGTGGTGAACA		

Table 1. Primer sequences of the gene-specific primers for different targets in the ICEPlex Fungal Assay

Results

Simultaneous detection of fungal targets in one reaction





Figure 3. Simultaneous detection of single- and multi-copy targets in the ICEPlex Fungal Assay.



The ICEPlex Fungal Assay can simultaneously detect and differentiate six major fungal pathogens. Each fungal species is detected by two amplicons: one is from single-copy region and the other from multi-copy region on chromosomal or mitochondrial genome (Fig.3). Based on calibration control copy number, assay can perform approximate quantification.

Embedded calibration, extraction and sensitivity controls




The ICEPlex Fungal Assay contains a sensitivity control (positive control for the presence of primers in the PCR reaction), an extraction control (*S. pombe*), and an optional Internal Control (data not shown here). The sensitivity control can be used either alone or in each PCR reaction.

Sensitivity and dynamic range of the ICEPlex Fungal Assay



Specificity of the ICEPlex Fungal Assay



Sample preparation method for whole blood

- Lysis**
 - 3 ml whole blood + 9 ml water and invert tube 5-10 times
 - 4,000 rpm for 2 min (tabletop)
 - Resuspend pellet in 9 ml water and invert tube 5-10 times
 - 4,000 rpm for 2 min (tabletop)
 - Resuspend pellet in 1 ml Protease Solution (30 μl Proteinase K in RBC Lysis Buffer), 15 min at 60 °C
 - 14,000 rpm for 2 min (microcentrifuge)
 - Resuspend pellet in 0.2 ml Lyticase Solution (20U Lyticase in 50mM Tris-Cl, pH 7.6, 1 mM EDTA, 0.2% β-mercaptoethanol), 30 min at 37 °C
 - No centrifugation
- Extraction**
 - Add 0.3ml Bead-Beating Buffer (50mM Tris-HCl, pH 7.6, 25 mM EDTA, 0.2% β-mercaptoethanol) for bead-beating
 - 4 m/sec for 60 sec on FastPrep-24
 - Harvest the lysate and rinse the beads with 0.6 5ml buffer for a final 1ml sample volume
 - DNA isolation on easyMAG (Specific B/Whole Blood Protocol), 50 μl elution

Conclusions

- The ICEPlex Fungal Assay consists of yeast-specific primers and primers for all controls, including a calibration control, a sensitivity control, an internal control, and an extraction control.
- No cross-reactivity was detected with gram-negative bacteria, gram-positive bacteria, other fungal genomic DNA mixes or human genomic DNA.
- No major multiplex primer interactions were observed. The PCR efficiencies of each fungal target primer set were comparable between single and multiplex formats (data not shown).
- Primers and PCR chemistry could detect single copies of DNA in PCR reaction. For most targets, assay showed five-dynamic range with PCR efficiencies ranging from 70-110%.
- Sample preparation method indicated effective removal of red and white blood cells. Further, fungal cell lysis method indicated a lysis efficiency of 90-100% and DNA extraction efficiency of 40-60% (data not shown).
- Study demonstrated the capability of the ICEPlex system to simultaneously detect major yeast pathogens in a single-tube multiplex format.

The ICEPlex system and ICEPlex Fungal Assay have not been approved for IVD by the FDA. The presented information is for demonstration purpose only.