

Application Brief

NRAS/BRAF Point Mutation Analysis Panel

Unique All-In-One-Well Assay with High Specificity and Sensitivity


INTRODUCTION

The *NRAS* and *BRAF* genes are proto-oncogenes involved in normal cellular function of controlling cellular growth. However, *NRAS* and *BRAF* point mutations can lead to tumor development and proliferation and have been implicated in variety of human tumors and particularly in melanoma. *NRAS* and *BRAF* somatic mutations are found in 9-29% and 53-66% of melanomas, respectively. These mutations result in continuous cell growth and are potential targets for therapy. Information on the mutation status can help clinicians choose appropriate treatment. Here we have developed a single-well multiplex *NRAS/BRAF* Point Mutation Analysis Panel on the ICEPlex® system that can detect and discriminate 12 *NRAS* and 4 *BRAF* clinically important mutations as shown below:

| CDS Mutation | Amino Acid |
|---------------------|------------|
| NRAS c.35 G>A | G12D |
| NRAS c.34 G>A | G12S |
| NRAS c.38 G>C | G13A |
| NRAS c.38 G>A | G13D |
| NRAS c.37 G>C | G13R |
| NRAS c.38 G>T | G13V |
| NRAS c.37 G>T | G13C |
| NRAS c.183 A>C (H1) | Q61H |

| CDS Mutation | Amino Acid |
|---------------------------|------------|
| NRAS c.182 A>G (R1) | Q61R |
| NRAS c.182_183 AA>GG (R2) | Q61R |
| NRAS c.182 A>T | Q61L |
| NRAS c.181 C>A | Q61K |
| BRAF c.1799_1800 TG>AT | V600D |
| BRAF c.1799 T>A | V600E |
| c.1799_1800 TG>AA | V600E |
| c.1798_1799 GT>AA | V600K |

MULTIPLE ANSWERS IN LESS THAN 4 HOURS



SUMMARY

- Delivers multiple answers - 16 *NRAS/BRAF* mutations in a single PCR reaction.
- Provides high quality by simultaneous detection of built-in controls such as: a set of calibration controls or size standards, DNA quality and extraction controls.
- Requires minimum DNA input and addresses specimen size issue.
- Expedites sample turn-around time to less than 4 hours.

METHOD HIGHLIGHTS

- Primers were designed using PrimeraDx's unique strategy that can selectively amplify *NRAS*/*BRAF* point mutations. All primers were analyzed in silico for primer-primer interaction and cross-reactivity. One of the primers in each primer set was labeled with either FAM- or TYE- dye.
- PCR amplification conditions were optimized using proprietary PCR chemistry on the ICEPlex system.
- Multiplex PCR reactions were subjected to thermocycling on a standard 96-well PCR plate on the ICEPlex system.
- The fluorescently labeled amplicons for the different *NRAS*/*BRAF* mutations were injected, separated and detected in the capillary electrophoresis module of the ICEPlex system.
- ICEPlex system software plotted the fluorescent signals for different amplicons, generated amplification curves for all targets and controls, and calculated cycle thresholds (Cts).

TYPICAL DATA

As shown below, we were able to detect and discriminate 13 important mutations in *NRAS* and 3 important mutations in *BRAF* genes.




Figure 1 Representative amplification curves for 3 *NRAS*/*BRAF* targets on the ICEPlex system.

| | | G12D | G12S | G13A | G13C | G13D | G13R | G13V | Q61H1 | Q61H2 | Q61K | Q61L | Q61R1 | Q61R2 | V600D | V600E | V600K |
|----------|--------|------|------|------|------|------|------|----------|-------|-------|------|----------|----------|-------|-------|-------|-------|
| Sample 1 | Ct | | | | | | | | | | | 26.2 | | | | | |
| | Result | No | No | No | No | Detected | No | No | No | No | No |
| Sample 2 | Ct | | | | | | | 23.5 | | | | | | | | | |
| | Result | No | No | No | No | No | No | Detected | No | No | No | No | No | No | No | No | No |
| Sample 3 | Ct | | | | | | | | | | | | 27.2 | | | | |
| | Result | No | No | No | No | No | Detected | No | No | No | No |

Figure 2. Representative results for three samples.

FOR MORE INFORMATION

For a list of publications and to find out more about how PrimeraDx can help your lab, please contact us at 508.618.2300 or visit www.primeradx.com.

The ICEPlex system and ICEPlex KRAS/BRAF Assay are for Research Use Only and have not been approved for in vitro diagnostic use by the FDA. The presented information is for demonstration purposes only.